

THE DENTAL ADVISOR

"Improving Patient Care Through Research & Education"

EDITORS

John W. Farah, D.D.S., Ph.D.
John M. Powers, Ph.D.

ADVISORY BOARD

John O. Burgess, D.D.S., M.S.
Gerard Kugel, D.M.D., M.S.
Karl F. Leinfelder, D.D.S., M.S.
John A. Molinari, Ph.D.
Frederick A. Rueggeberg, D.D.S., M.S.
Edward J. Swift, Jr., D.M.D., M.S.
Charles W. Wakefield, D.D.S.

DIRECTOR OF RESEARCH

Ron Yapp, M.S.

CREATIVE DIRECTOR

Lynn St. Pierre

EDITORIAL STAFF

Stephanie D. Farquhar
Annette M. Frederick
Tricia G. Price

ADMINISTRATIVE STAFF

Kathryn L. Jackson
Samantha Baker
Darlene E. Meyer

PUBLISHER

Dental Consultants, Inc.

THE DENTAL ADVISOR

3110 West Liberty
Ann Arbor, MI 48103

Toll free: 800.347.1330

Local: 734.665.2020

E-mail: info@dentaladvisor.com

Web Site: www.dentaladvisor.com

NUMBER 4 • JANUARY 2006

Effects of Preheating Resin Composite on Flowability

James C. Broome, D.D.S., M.S., University of Alabama at Birmingham School of Dentistry

Purpose – Preheating of resin composite compules has been suggested as a means of improving flow characteristics. The purpose of this study was to quantify in vitro the effect of preheating of resin composite compules on flowability.

Methods – A screening test was devised to assess the flow of six resin composites at room temperature (24° C) and after heating to 55° C and to 60° C in a commercial heating unit (*Calset/AdDent, Inc.*). Composite compules were placed into a steel positioning stand for heating and testing. A load was applied to the plunger of the compule by a steel rod mounted in a chuck on a servo-hydraulic testing machine (MTS 858, MTS Systems Corp.) at a crosshead speed of 1.0 mm/sec. The load (N) after displacement of 6 mm was recorded as a measure of flowability. Data were analyzed by analysis of variance and means were compared by a Fisher's PLSD test at the 0.05 level of significance.

Results – Preliminary tests indicated that 1.0 mm/sec was the approximate rate that a load would be applied a compule with a conventional injection gun. Heating in a commercial unit produced a statistically significant improvement in flow ranging from 5 to 76% for the resin composites tested.

Table 1. Load (N) at a displacement of 6 mm for each material at each temperature*.

	<i>Esthet-X</i> DENTSPLY Caulk <i>Microhybrid</i>	<i>Filtek Supreme</i> 3M ESPE <i>Nanofilled</i>	<i>Heliomolar</i> Ivoclar Vivadent <i>Microfilled</i>	<i>Point 4</i> SDS/Kerr <i>Microhybrid</i>	<i>Simile</i> Pentron Clinical Technologies, LLC <i>Nanofilled</i>	<i>Filtek Z250</i> 3M ESPE <i>Microfine</i>
24° C	215 (11)	94.3 (8.5)	81.9 ^b (9.2)	98.7 (7.3)	105 (19)	191 (21)
55° C	108 ^a (13)	51.1 (5.7)	77.6 ^b (6.3)	40.4 ^c (3.3)	68.7 (6.1)	57 ^d (13)
60° C	96 ^a (14)	39.5 (4.4)	63.6 (6.3)	36.4 ^c (3.2)	47.4 (2.3)	45.0 ^d (4.3)

* Means with standard deviations in parentheses (n=5). Means with the same superscript were not significantly different (p>0.05).

Table 2. Increase in flow (%) for each resin composite from baseline at 55° and at 60° C.

	<i>Esthet-X</i>	<i>Filtek Supreme</i>	<i>Heliomolar</i>	<i>Point 4</i>	<i>Simile</i>	<i>Filtek Z250</i>
55° C	50	46	5	59	34	70
60° C	55	58	22	63	55	76

Conclusion – Preheating resin composite compules at 55° and 60° C provided an effective means of increasing flowability of the composites tested.

Acknowledgment – This project was supported in part by AdDent through a grant to the University of Alabama at Birmingham School of Dentistry.