

### **COMPLUES THAT FIT THE COMPEX DISPENSER**

The Complex was designed to fit most compules on the market with a diameter of 6.32 mm to 6.68mm.

<b>No.</b>	<b>Manufacturer</b>	<b>Composite Name</b>
1	3M ESPE	Filtek Supreme
2	Clinician's Choice	Evanisce
3	Cosmedent	Renamel Nano
4	Dentsply	TPH Spectra
5	DMG America	Ecosite Bulk Fill
6	GC	G-aenial Sculpt
7	Ivoclar	Tetric EvoCeram
8	Kerr	Harmonize
9	Kulzer	Venus
10	Kuraray	AP-X
11	Septodont	N'Durance
12	Shofu	Beautifil
13	Tokuyama	Estelite
14	Ultradent	Amelogen plus
15	VOCO	Admira Fusion
16	VOCO	GrandioSO
17	VOCO	X-tra Fil
18	VOCO	VisCalor
19	VOCO	VisCalor bulk



**After continuous use of larger diameter compules, check to ensure that smaller diameter compules have a secure fit in the compule holder prior to use.**