

# Compex HD™

THERMAL ASSISTED  
COMPOSITE DISPENSER



43 Miry Brook Road • Danbury, CT 06810 • (203) 778-0200

**AdDent**  
INCORPORATED

Advancing Dental  
Excellence

ISO 9001 CERTIFIED & 13485 COMPLIANT | ALL PRODUCTS ARE FDA CLEARED

Made in the USA | [www.addent.com](http://www.addent.com)

## COMPEX HD™

The Compex HD™ Uses the Latest PTC Technology to Rapidly Warm Heavily Filled Compules

AdDent, Inc., the first to market composite warmers, has re-envisioned composite delivery with the Compex HD™ Heated Dispenser. This innovative device uses PTC technology to rapidly warm even heavily filled compules.

- Able to quickly heat composite material to 155°F
- Warming composite significantly lowers the viscosity of the material, resulting in better adaptation, reduced voids, gaps, and microleakage
- Reduces cure time and increases the depth of cure
- Increases micro-hardness and reduces shrinkage stress
- Materials remain highly sculptable, non-sticky, and easily shaped
- Improves shade stability

### FEATURES

- Light weight handheld rechargeable device.
- Heats a broad range of compules.
- Heats compules to 155°F (68°C) in 40 seconds or less.
- Heats up to 100 compules without recharging.
- Easy to load and unload with patented ejector mechanism.
- Silicone sleeve and barrier bag.
- Short recharge time.

### BENEFITS

- Comfortable all in one heater and dispenser – saves space, easily moved room to room.
- Heats most composite compules currently on the market.
- Delivers all composite without temperature loss.
- Charge no more than once per week under normal use.
- Fast compule insertion and removal saves time.
- Minimizes cross contamination by dentist or assistant.
- Full charge time is only 4 ½ hours.

## COMPEX HD™ WILL ENHANCE PLACEMENT FOR ANY PROTOCOL

The Compex HD™ is extremely versatile and fits most brands of composite compules. This device enhances posterior bulk fills and traditional incremental layering techniques. Delivering all the benefits of heated composite, directly to the restoration by maintaining the optimal temperature during the entire protocol.

## STUDIES PROVE THE BENEFITS OF HEATING COMPOSITE

The heating technology of Compex HD™ lowers the viscosity of all composite materials, allowing it to better flow and adapt to preparation walls. As a result, voids and micro-leakage are reduced, which then reduces the risk of secondary caries and restorative failures. Heated composite material also improves the physical and handling properties while improving color stability.

References available at [www.addent.com](http://www.addent.com)

Using the Compex HD™ for all your composite protocols, allows you to experience procedure efficiency. Placing all the benefits of heated composite in your hand

## COMPEX HD™ HEATED DISPENSER:

Includes: Compex HD™ (1) dispenser, (10) silicone sleeves, (25) barrier bags, and (1) medical grade charge with USB cord

Follow us on social media:



[info@addent.com](mailto:info@addent.com)



SCAN ME

**AdDent**  
INCORPORATED

Advancing Dental  
Excellence

ISO 9001 CERTIFIED & 13485 COMPLIANT | ALL PRODUCTS ARE FDA CLEARED

Made in the USA | [www.addent.com](http://www.addent.com)